

Poly(methyl methacrylate) (PMMA) XPS Reference Core Level and Energy Loss Spectra

Pierre Louette, Frederic Bodino, and Jean-Jacques Pireaux

Citation: [Surface Science Spectra](#) **12**, 69 (2005); doi: 10.1116/11.20050914

View online: <http://dx.doi.org/10.1116/11.20050914>

View Table of Contents: <http://scitation.aip.org/content/avs/journal/sss/12/1?ver=pdfcov>

Published by the AVS: Science & Technology of Materials, Interfaces, and Processing

Articles you may be interested in

[Poly\(acroleine\) XPS Reference Core Level and Energy Loss Spectra](#)

Surf. Sci. Spectra **12**, 154 (2005); 10.1116/11.20051107

[Poly\(butyl methacrylate\) \(PBMA\) XPS Reference Core Level and Energy Loss Spectra](#)

Surf. Sci. Spectra **12**, 139 (2005); 10.1116/11.20051104

[Poly\(methyl hydrogen siloxane\) \(PMHS\) XPS Reference Core Level and Energy Loss Spectra](#)


Surf. Sci. Spectra **12**, 133 (2005); 10.1116/11.20051103

[Poly\(styrene\) XPS Reference Core Level and Energy Loss Spectra](#)




Surf. Sci. Spectra **12**, 96 (2005); 10.1116/11.20050920

[Poly\(ethyl methacrylate\) \(PEMA\) XPS Reference Core Level and Energy Loss Spectra](#)

Surf. Sci. Spectra **12**, 44 (2005); 10.1116/11.20050909



Instruments for Advanced Science

<p>Contact Hiden Analytical for further details: W www.HidenAnalytical.com E info@hiden.co.uk CLICK TO VIEW our product catalogue</p>	 <p>Gas Analysis</p> <ul style="list-style-type: none"> › dynamic measurement of reaction gas streams › catalysis and thermal analysis › molecular beam studies › dissolved species probes › fermentation, environmental and ecological studies 	 <p>Surface Science</p> <ul style="list-style-type: none"> › UHV-TPD › SIMS › end point detection in ion beam etch › elemental imaging - surface mapping 	 <p>Plasma Diagnostics</p> <ul style="list-style-type: none"> › plasma source characterization › etch and deposition process reaction › kinetic studies › analysis of neutral and radical species 	 <p>Vacuum Analysis</p> <ul style="list-style-type: none"> › partial pressure measurement and control of process gases › reactive sputter process control › vacuum diagnostics › vacuum coating process monitoring
---	--	--	--	--

Poly(methyl methacrylate) (PMMA) XPS Reference Core Level and Energy Loss Spectra

Pierre Louette,^{a)} Frederic Bodino, and Jean-Jacques Pireaux
University of Namur - LISE, 61, Rue de Bruxelles, B-5000 Namur, Belgium

(Received 23 September 2005; accepted 6 June 2006; published 30 January 2007)

XPS measurements of poly(methyl methacrylate) recorded with a SSX-100 spectrometer in standardized experimental conditions are presented: survey scan, high resolution core level spectra as well as the energy loss regions of carbon and oxygen peaks are analyzed. This is part of a contract work aiming to record spectra in the very same conditions of some 40 different polymers. © 2006 American Vacuum Society. [DOI: 10.1116/11.20050914]

Keywords: *x-ray photoelectron spectroscopy; XPS; surface; polymer; poly(methyl methacrylate); PMMA*

PACS: 79.60.Fr, 82.80.Pv, 79.20.Uv, 61.41.+e, 01.30.Kj

Accession # 00805

Technique: XPS

Host Material: poly(methyl methacrylate) (PMMA)

Instrument: Surface Science Instruments SSX-100

Major Elements in Spectrum: C, O

Minor Elements in Spectrum: none

Printed Spectra: 5

Spectra in Electronic Record: 5

Spectral Category: comparison

SPECIMEN DESCRIPTION

Host Material: poly(methyl methacrylate) (PMMA)

CAS Registry #: 9011-14-7


Host Material Characteristics: homogeneous; solid; amorphous; dielectric; polymer

Chemical Name: poly(methyl methacrylate)

Source: Scientific Polymer Products, Inc.

Host Composition: (C₅H₈O₂)_n

Form: beads



History & Significance: This study is a part of a reference spectra database of polymers, including survey and core level spectra, but also energy loss spectra of the main elements.

As Received Condition: not specified

Analyzed Region: same as host material

Ex Situ Preparation/Mounting: not specified

In Situ Preparation: none

Pre-Analysis Beam Exposure: The analyzed region was exposed to x-rays for a very short time, around 2 min for sample position adjustment prior to measurements.

Charge Control: use of a metal screen and a flood gun (6 eV)

Temp. During Analysis: 300 K

Pressure During Analysis: <1.2×10⁻⁶ Pa

INSTRUMENT DESCRIPTION

Manufacturer and Model: Surface Science Instruments SSX-100

Analyzer Type: spherical sector

Detector: position sensitive detector with microchannel plate

Number of Detector Elements: 128

INSTRUMENT PARAMETERS COMMON TO ALL SPECTRA

■ Spectrometer

Analyzer Mode: constant pass energy

Throughput ($T = E^N$): $N = \text{See comment below}$

Throughput Comment: $T = E^N$, $N = 0.7$

Excitation Source Window: 1.5 μm Al foil

Excitation Source: Al K_α monochromatic

Source Energy: 1486.6 eV

Source Strength: 130 W

Source Beam Size: 0.6 mm × 0.6 mm

Signal Mode: not specified

■ Geometry

Incident Angle: 57.6°

Source to Analyzer Angle: 70.8°

Emission Angle: 14.7°

Specimen Azimuthal Angle: 75.5°

Acceptance Angle from Analyzer Axis: 0°

Analyzer Angular Acceptance Width: 30° × 30°

DATA ANALYSIS METHOD

Energy Scale Correction: To compensate for charging effects, we adjusted the largest C1s component to 285.00 eV (Ref. 1).

Recommended Energy-Scale Shift: +9.41 eV

Peak Shape and Background Method: A least square fitting routine with mixed Gaussian/Lorentzian for the components and a linear background was used.

Quantitation Method: Scofield factors corrected for energy dependence were used.

^{a)}Author to whom correspondence should be addressed.

ACKNOWLEDGMENTS


This study is a part of the EU-BCR contract “XPS Spectral Intensity Data Bank”. We thank the NPL for authorizing us to publish these spectra.

REFERENCES

1. G. Beamson and D. Briggs, in *The Scienta ESCA 300 Database* (Wiley, Chichester, 1992).
2. C. J. Powell, J. Electron. Spectrosc. Relat. Phenom. **47**, 197 (1988).

SPECTRAL FEATURES TABLE							
Spectrum ID #	Element/Transition	Peak Energy (eV)	Peak Width FWHM (eV)	Peak Area (eV-cts/s)	Sensitivity Factor	Concentration (at. %)	Peak Assignment
00805-02	C 1s	285.00	1.30	3169	1.00	31.7	1 in Diagram below
00805-02	C 1s	285.84	1.30	1471	1.00	14.7	2 in Diagram below
00805-02	C 1s	286.81	1.30	1303	1.00	13.1	3 in Diagram below
00805-02	C 1s	289.00	1.21	1313	1.00	13.2	4 in Diagram below
00805-03	O 1s	532.28	1.40	3322	2.49	13.4	5 in Diagram below
00805-03	O 1s	533.79	1.57	3450	2.49	13.9	6 in Diagram below


Comment to Spectral Features Table:




ANALYZER CALIBRATION TABLE							
Spectrum ID #	Element/Transition	Peak Energy (eV)	Peak Width FWHM (eV)	Peak Area (eV-cts/s)	Sensitivity Factor	Concentration (at. %)	Peak Assignment
...	Cu 2p _{3/2}	932.34	1.19	202906	9.748
...	Cu 3p _{3/2}	74.78	2.36	289045	2.774

GUIDE TO FIGURES					
Spectrum (Accession) #	Spectral Region	Voltage Shift*	Multiplier	Baseline	Comment #
805-1	Survey	−9.41	1	0	
805-2	C 1s	−9.41	1	0	
805-3	O 1s	−9.41	1	0	
805-4	C 1s + losses	−9.41	1	0	
805-5	O 1s + losses	−9.41	1	0	

* Voltage shift of the archived (as-measured) spectrum relative to the printed figure. The figure reflects the recommended energy scale correction due to a calibration correction, sample charging, flood gun, or other phenomenon.




Accession #	00805-01
Host Material	poly(methyl methacrylate) (PMMA)
Technique	XPS
Spectral Region	survey
Instrument	Surface Science Instruments SSX-100
Excitation Source	Al K_{α} monochromatic
Source Energy	1486.6 eV
Source Strength	130 W
Source Size	0.6 mm \times 0.6 mm
Analyzer Type	spherical sector
Incident Angle	57.6°
Emission Angle	14.7°
Analyzer Pass Energy	106.8 eV
Analyzer Resolution	1.17 eV
Total Signal Accumulation Time	960 s
Total Elapsed Time	not specified
Number of Scans	2
Effective Detector Width	12.96 eV




■ **Accession #:** 00805-02
 ■ **Host Material:** poly(methyl methacrylate) (PMMA)
 ■ **Technique:** XPS
 ■ **Spectral Region:** C 1s

Instrument: Surface Science Instruments SSX-100
 Excitation Source: Al K_{α} monochromatic
 Source Energy: 1486.6 eV
 Source Strength: 130 W
 Source Size: 0.6 mm \times 0.6 mm
 Incident Angle: 57.6°
 Analyzer Type: spherical sector
 Analyzer Pass Energy: 29.97 eV
 Analyzer Resolution: 0.76 eV
 Emission Angle: 14.7°
 Total Signal Accumulation Time: 1200 s
 Total Elapsed Time: not specified
 Number of Scans: 10
 Effective Detector Width: 3.341 eV




■ **Accession #:** 00805-03
 ■ **Host Material:** poly(methyl methacrylate) (PMMA)
 ■ **Technique:** XPS
 ■ **Spectral Region:** O 1s

Instrument: Surface Science Instruments SSX-100
 Excitation Source: Al K_{α} monochromatic
 Source Energy: 1486.6 eV
 Source Strength: 130 W
 Source Size: 0.6 mm \times 0.6 mm
 Incident Angle: 57.6°
 Analyzer Type: spherical sector
 Analyzer Pass Energy: 29.97 eV
 Analyzer Resolution: 0.76 eV
 Emission Angle: 14.7°
 Total Signal Accumulation Time: 1200 s
 Total Elapsed Time: not specified
 Number of Scans: 10
 Effective Detector Width: 3.341 eV



Accession #: 00805-04
Host Material: poly(methyl methacrylate) (PMMA)
Technique: XPS
Spectral Region: C 1s energy losses
 Instrument: Surface Science Instruments SSX-100
 Excitation Source: Al K_{α} monochromatic
 Source Energy: 1486.6 eV
 Source Strength: 130 W
 Source Size: 0.6 mm \times 0.6 mm
 Incident Angle: 57.6°
 Analyzer Type: spherical sector
 Analyzer Pass Energy: 106.8 eV
 Analyzer Resolution: 1.17 eV
 Emission Angle: 14.7°
 Total Signal Accumulation Time: 900 s
 Total Elapsed Time: not specified
 Number of Scans: 5
 Effective Detector Width: 12.956 eV



Accession #: 00805-05
Host Material: poly(methyl methacrylate) (PMMA)
Technique: XPS
Spectral Region: O 1s energy losses
 Instrument: Surface Science Instruments SSX-100
 Excitation Source: Al K_{α} monochromatic
 Source Energy: 1486.6 eV
 Source Strength: 130 W
 Source Size: 0.6 mm \times 0.6 mm
 Incident Angle: 57.6°
 Analyzer Type: spherical sector
 Analyzer Pass Energy: 106.8 eV
 Analyzer Resolution: 1.17 eV
 Emission Angle: 14.7°
 Total Signal Accumulation Time: 900 s
 Total Elapsed Time: not specified
 Number of Scans: 5
 Effective Detector Width: 12.956 eV